

# FreeWheel ADAPTOR for *Folding Wheelchair*

Individuals with folding wheelchairs may now experience the independence, freedom and added safety the FreeWheel provides.


Using the horizontal side tubing, the Adaptor securely clamps on to a folding wheelchair. Once set up for the individual chair, it comes on and off in seconds, making it easy to use.

The Adaptor and FreeWheel may be used on wheelchairs that have 2 or 1 flip up footrest types. Maximum weight of end user, 250 lbs. or less.

**Freedom. Independence. Mobility.**

## SPECIFICATIONS

- ◆ Need 2.5" of round, horizontal side tubing on each side of wheelchair to attach Adaptor. May be angled, as long as the sides are parallel.
- ◆ Round horizontal tubing may be 1" or 1-1/8" in diameter. Custom tube sleeves available for 3/4" and 7/8" tubing.
- ◆ Width between horizontal tubing 13" to 19". Customization options available for small width (pediatric wheelchairs)
- ◆ Side tubing height: 2 footrests 6' to 16". Single footrest 10" to 16" with footrest height max of 4" (custom frame options may exist).
- ◆ FreeWheel clamps on to "HITCH" using custom Thick Flat shim and Short Rear End.
- ◆ Easily remove Adaptor by unscrewing knobs on side tubing. Knobs may be positioned on the top or underneath side tubing.
- ◆ Maximum personal body weight for ADAPTOR use is 250 lbs.
- ◆ Weight of ADAPTOR is 3.5 lbs.



For more information please contact us at:

Email: [info@gofreewheel.com](mailto:info@gofreewheel.com)

Phone: 208-571-2051

[www.gofreewheel.com](http://www.gofreewheel.com)